

The *Streptococcus pyogenes* MLST scheme uses internal fragments of seven housekeeping genes amplified by PCR using the following primer pairs (described in [Enright MC, Spratt BG, Kalia A, Cross JH and Bessen DE. 2001. Multilocus sequence typing of *Streptococcus pyogenes* and the relationships between *emm* type and clone. Infection and Immunity 69:2416-2427](#)):

Genes and Function	Sequences (5' to 3')	Size (bp) of the amplicon used for assigning alleles
glucose kinase		
<i>gki-up</i>	GGCATTGGAATGGGATCACC	498
<i>gki-dn</i>	TCTCCTGCTGCTGACAC	
glutamine transporter protein		
<i>gtr-up</i>	GAGGTTGTGGTGATTATTGG	450
<i>gtr-dn</i>	GCAAAGCCCATTCATGAGTC	
glutamate racemase		
<i>murI-up</i>	TGCTGACTCAAATGTTAAAATGATTG	438
<i>murI-dn</i>	GATGATAATTCACCGTTAATGTCAAATAG	
DNA mismatch repair protein		
<i>mutS-up</i>	GAAGAGTCATCTAGTTTAGAATACGAT	405
<i>mutS-dn</i>	AGAGAGTTGTCACTTGCGCGTTTGATTGCT	
transketolase		
<i>recP-up</i>	GCAAATTCTGGACACCCAGG	459
<i>recP-dn</i>	CTTTCACAAGGATATGTTGCC	
xanthine phosphoribosyl transferase		
<i>xpt-up</i>	TTACTTGAAGAACGCATCTTA	450
<i>xpt-dn</i>	ATGAGGTCACCTCAATGCC	
acetyl-CoA acetyltransferase		
<i>yqiL-up</i>	TGCAACAGTATGGACTGACCAGAGAACAAGATGC	434
<i>yqiL-dn</i>	CAAGGTCTCGTGAAACCGCTAAAGCCTGAG	

PCR conditions

The PCR reactions are performed in volumes of 50 μ L, with an initial denaturation at 95°C for five min, followed by 28 cycles of 95°C for 1 min, 55°C for 1 min and 72°C for 1 min. The amplified DNA fragments are purified either by precipitation with polyethylene glycol or using a commercial PCR purification kit. The sequence of each fragment is obtained on both strands using the same primers as those in the initial PCR amplifications.

As the same primers are used for amplification and sequencing, it is important that only a single DNA fragment is amplified in the initial PCR. This may involve some optimisation of the annealing temperature and other PCR conditions in individual laboratories.